**TMDA/DMD/MDA/F/014.**  
**Rev #:00**

**THE UNITED REPUBLIC OF TANZANIA**

**MINISTRY OF HEALTH**

**TANZANIA MEDICINES AND MEDICAL DEVICES AUTHORITY**

**PUBLIC ASSESSMENT REPORT FOR WONDFO HIV SELF-TEST**

**Version number 2.0, 29/03/2024**

**TMDA Headquarters, Plot No. 56/1, Block E, Kisasa B Centre, Hombolo Road, P. O. Box 1253, Dodoma – Tanzania, Telephone: +255 (26) 2961989/2061990/+255 (22) 2450512/2450751/2452108, Email: [info@tmda.go.tz](mailto:info@tmda.go.tz), Website: [www.tmda.go.tz](http://www.tmda.go.tz)**

**Toll free: 0800110084**

## 1. Introduction

Wondfo HIV Self-Test is a class D in-vitro diagnostic device belonging to the microbiology specialty category. Wondfo HIV Self-Test is approved in Tanzania as a kit, for use by the general public.

### 1.1. Administrative Information

Registration number	TAN 23 MDR 0157
Brand Name (if relevant)	Wondfo HIV Self-Test
Common name	HIV Self-Test
Class of the device and rule applied	Class D according to Rule 1 of Classification for In Vitro Diagnostic Devices
GMDN code and term	65848 HIV1/HIV2 antibody IVD, kit, rapid ICT, self-testing
Name and complete address of the Market Authorization Holder	Guangzhou Wondfo Biotech Co., Ltd. No. 8 Lizhishan Road, Science City, Huangpu District, 510663 Guangzhou, P.R. China. Telephone: +86 18565352901, Fax: 86-20-32296062/32296063, E-mail: <a href="mailto:global@wondfo.com.cn">global@wondfo.com.cn</a> , Website: <a href="http://en.wondfo.com">en.wondfo.com</a> .
Name and address(es) of local responsible person (LRP).	Horizon Pharmacy Ltd, Faykat towers, Block 41, Morocco, Kinondoni, P.O.BOX 31405, Dar es salaam, Tanzania, Telephone: +255 715 268781, E-mail: <a href="mailto:horizonpharmacy22@gmail.com">horizonpharmacy22@gmail.com</a> ,

### 1.2. Assessment Procedure

The application for registration of Wondfo HIV Self-Test was submitted on 15/05/2023. The product underwent full procedure assessment. Assessment was completed in 03 rounds of evaluation. Wondfo HIV Self-Test was registered on 13/07/2023.

## 2. Technical information

### 2.1. Intended use

The intended use of Wondfo HIV Self-Test as declared by the manufacturer and approved by TMDA is screening for the detection of antibodies to HIV-1 and HIV-2 in blood to aid in the diagnosis of infection with HIV-1 and HIV- 2. Wondfo HIV Self-Test is approved for use in point of care by both trained personnel and laypersons.

### 2.2. Device details and features

Wondfo HIV Self-Test has been registered as a kit which consists of test cassette, buffer, alcohol wipe, sterile pad, plaster, and a single-use lancet.

Wondfo HIV Self-Test is an in vitro diagnostic device. It is used for screening of HIV-1 and HIV- 2 infection. Wondfo HIV Self-Test operates by immunochromatographic principle. The test out-put is qualitative.


The type of specimen used is whole blood and is collected by pricking a finger using a lancet provided in the kit to obtain a drop of blood.

#### Device description

The Wondfo HIV Self-Test adopts double antigen sandwich immune chromatography method. Once the specimen and the buffer are added into their respective wells, they migrate along the device by capillary action. When the levels of HIV-1 or HIV-2 antibodies are at or above the target cut-off level, HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41). The complex is captured by HIV recombinant antigen (gp41and gp36) immobilized in the test region (T). This produces a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV 1/2 antibody are zero or below the target cut- off level, there is no visible colored band in the test region (T). This indicates a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C), if the test has been performed properly.

#### Pictorial Presentation

## OPERATION PROCEDURE

### 2.3. Commercial presentation

There are 03 approved commercial presentations as follows: one test cassette in a pouch with desiccant. 1 or 20 or 100 pouches are placed in secondary packaging material.

additional contents include.

- 1 or 20 or 100 Buffers
- 1 or 20 or 100 Droppers
- 1 or 20 or 100 Blood lancets for single use
- 1 or 20 or 100 Cotton swabs
- 1 or 20 or 100 Alcohol prep pads
- 1 disposable bag (Only for 1 Test/kit)
- 1 Or 20 or 100 Instruction for use

### 2.4. Items required but not submitted,

- Timer
- Tissue
- Bandage

## 3. Storage instructions

### 3.1.1. Shelf-life

The approved shelf-life is 24 months.

### 3.1.2. Storage conditions

The recommended storage conditions are at 2-30 °C as declared by the manufacturer. Use the kit within one hour after opening.

### 3.1.3. Shipping conditions

The recommended shipping conditions are 2-30 °C as declared by the manufacturer.

## 4. Manufacturing site audit

The manufacturer of the device is,  
Guangzhou Wondfo Biotech Co., Ltd,  
No. 8 Lizhishan Road, Science City, Huangpu District,  
510663 Guangzhou, P.R. China.

Quality audit of the manufacturing facility was conducted through site visit on 12<sup>th</sup> – 13<sup>th</sup> June, 2019. The site was found to be compliant with ISO 13485 requirements.

## 5. Performance Evaluation

### 5.1. Analytical Performance

The analytical performance characteristics of the device was established through the following test parameters: accuracy, trueness, precision (repeatability and reproducibility), sensitivity, and specificity.

### 5.2. Clinical Performance

Clinical performance was conducted at; National Institute of Medical Research (NIMR), Mwanza Medical Research Centre, P.O. Box 11936, Isamilo St., Ilemela – Mwanza, Tel: +255 (0) 28 250 0019, [info@mitu.or.tz](mailto:info@mitu.or.tz),

The following parameters were tested sensitivity and specificity. Based on results of the performance studies, it was concluded that the test sensitivity and specificity is 100% and 99.96% for professional users. The test sensitivity and specificity are 95.8% and 99.6% for Self-Test users respectively. The studies further concluded that Wondfo HIV Self-Test is capable of consistently producing accurate and reliable test output.

## 6. Product label and instructions for use

The content of the primary and secondary labels is in line with TMDA labeling requirements in terms of content, layout and design. The label contains sufficient information for proper identification of the device and post marketing follow up of the product.


The instruction for use includes all the relevant information to ensure correct and safe use of the device by both trained personnel and lay persons.

6.1. Primary pack


6.2. Secondary pack

1. The box (outer packaging container) of Wondfo HIV Self-Test (20T/kit)

1. The box (outer packaging container) of Wondfo HIV Self-Test (100T/kit)

### 6.3 Instructions for use/Package insert.

Instructions for use can be accessed by clicking [here](#).

## 7. Risk – Benefit Analysis

On basis of the data submitted, the current state of knowledge and compliance of the manufacturer to ISO 13485, the benefit of the product outweighs the risks associated with its use when used in accordance with the manufacturer’s instruction. Wondfo HIV Self-Test was recommended for registration.

## 8. Post-approval updates

### 8.1. Variation applications

Reference number	Date submitted	Change requested	Recommendation	Granting date
NA				

**8.2. Feedback from pharmacovigilance, post marketing surveillance and enforcement activities**

Type of feedback	Impact	Response
No recorded Adverse Event		NA

**8.3. Re-registration applications**

NA

**CHANGE HISTORY**

Version number	Date	Description of update	Section(s) Modified	Approval date


# Wondfo HIV Self-Test


## Instruction for use

- You must follow the test procedure carefully to get an accurate result.
- Sit in a clean, well-lit area and ensure you have all contents before beginning the test.
- For single-use only. Do not open foil pouch containing cassette until ready to test.


## Contents

## How to use the test kit (for fingerstick whole blood use)

## How to read the test

### DISPOSE

Throw away all contents in the normal trash or dispose according to local infectious waste disposal laws and regulations.

## Questions and answers

### 1. What is HIV?

The human immunodeficiency virus (HIV) targets cells of the immune system, called CD4 cells, which help the body respond to infection. Within the CD4 cell, HIV replicates and in turn, damages and destroys the cell. Without effective treatment of a combination of antiretroviral (ARV) drugs, the immune system will become weakened to the point that it can no longer fight infection and disease.

### 2. How is HIV transmitted?

HIV is found in certain bodily fluids of people living with HIV, including blood, semen, vaginal fluids, rectal fluids and breastmilk. HIV can be transmitted by:

- Unprotected vaginal or anal sex, and, in very rare cases, through oral sex with a person living with HIV.
- Blood transfusion of contaminated blood.
- Sharing of needles, syringes, other injecting equipment, surgical equipment or other sharp instruments.
- From a mother living with HIV to her infant during pregnancy, childbirth or breastfeeding.

If a person living with HIV is on ART, which effectively suppresses HIV in the body, their chance of transmitting HIV to another person is greatly reduced.

### 3. When do I need to test myself?

Most widely-used HIV diagnostic tests detect antibodies produced by the person as part of their immune response to fight HIV. In most cases, people develop antibodies to HIV within 28 days of infection. During this time, people experience the so-called "window" period. Those reporting potential HIV exposure in the preceding 3 months following a negative self-test result need test again 3 months after the risk event to be sure.

### 4. What should I do if I get a positive result?

You need to follow up with a health care to get additional testing to confirm the result. At that time your local clinic and the doctor will suggest the next steps that need to be taken.

## Product information

### CATALOG NO.

W006P0058      W006P0059      W006P0060

### INTENDED USE

The Wondfo HIV Self-Test is a single-use *in vitro* diagnostic self-test for fingerstick whole blood detection of HIV-1/2. It is intended to be used by people who believe they may have been exposed to HIV or may have been at risk within the last 3 months. It is intended to be used as self-test and/or by medical professionals.

For *in vitro* diagnostic use only.

### PRINCIPLES OF THE PROCEDURE

Wondfo HIV Self-Test adopts double antigen sandwich immunochromatography method. Once the specimen and the buffer are added into the respective wells, they will migrate along the device by capillary action. The HIV-1 and/or HIV-2 antibodies bind to the colloidal gold-HIV antigen (gp36/41), and the complex is then captured by the HIV recombinant antigen (gp41 and gp36) immobilized in the test region (T). When the levels of HIV-1 or HIV-2 antibodies are at or above the limit of detection (LOD) of the assay, it will produce a visible colored band in the test region (T) and indicates a positive result. When the levels of HIV-1 or HIV-2 antibodies are zero or below the LOD, there is no visible colored band in the test region (T) indicating a negative result. To serve as an internal procedure control, a colored line will appear at the control region (C).

### WARNINGS AND CAUTIONS

**Do not** use if the test kit beyond expiration date.  
**Do not** use if the pouch is punctured or improperly sealed.  
**Do not** use for self-testing if you are under 12 years old.  
**Do not** use for self-testing if you have a bleeding disorder.  
**Do not** use for self-testing if you are already diagnosed as HIV positive.  
**Do not** open the pouch until you are ready to perform the test.  
 Wash your hands and ensure that they are clean and dry before starting test.  
 Adequate lighting is required to read the test results.

### KIT CONTENT

There are 3 configurations of the test kits, 1 test/kit, 20 tests/kit and 100 tests/kit. The kit components are provided as below:

Components		Catalog No.	W006P0058	W006P0059	W006P0060
Test cassette pouch	Test cassette (pcs)		1	1x20	1x100
	Desiccant (pcs)		1	1x20	1x100
	Dropper (pcs)		1	1x20	1x100
	Buffer (vial)		1	1x20	1x100
Accessories	IFU (pcs)		1	1x20	1x100
	Blood Lancet for Single Use (pcs)		1	1x20	1x100
	Alcohol prep pad (pcs)		1	1x20	1x100
	Cotton swab (pcs)		1	1x20	1x100
	Disposal bag (pcs)		1	/	/

One test strip includes: Gold conjugate (HIV gp41/gp36 fusion recombinant antigen-gold colloid and rabbit IgG polyclonal antibody-gold colloid), Test line (HIV gp41 recombinant antigen and HIV gp36 recombinant antigen) and Control line (Goat anti rabbit IgG polyclonal antibody).

### STORAGE AND STABILITY

- The test kit can be stored at 2-30 °C for 24 months.
- Use the test cassette within 1 hour after opening the pouch.
- Keep away from sunlight, moisture and heat.
- Use the kit at 10-30 °C.

### LIMITATIONS OF THE PROCEDURE

- The test is designed for detecting human fingerstick whole blood.
- The test is limited to the qualitative detection of HIV-1 and HIV-2 antibodies.
- The assay procedure and result interpretation must be followed closely when testing. For optimal test performance, proper specimen collection is critical. Failure to follow the procedure may lead to inaccurate test results.
  - Patients exposed to HIV less than 3 months.
  - Patients under HIV treatment (Anti retrovirus drugs).
  - If the quantity of antibodies for HIV present in the specimen is below the detection limit of the assay.
- False-negative results can occur in the following conditions:
  - Patients have participated in a HIV vaccine clinical trial.
- False-positive results can occur in the following conditions:
  - The correct specimen has been used.
  - The specimen has been applied correctly.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but instead should be determined by a healthcare provider in conjunction with clinical findings and the results from other laboratory tests and evaluations. Results from the Wondfo HIV Self-Test should not be used as the sole basis for diagnosis.

### PERFORMANCE CHARACTERISTICS

In the clinical study, 900 participants whose HIV status were unknown were given the Wondfo HIV Self-Test to test. The results were compared to the 4<sup>th</sup> generation laboratory test. The laboratory results shown that a total of 72 participants were HIV positive and 785 participants were HIV negative. 42 participants were excluded from the performance analysis. The comparison of results were as follows:

- 95.8% of participants (69/72) correctly reported the result as positive. This means that 3 participants infected with HIV reported negative result. This is called a false negative.
- 99.6% of participants (782/785) correctly report the result as negative. This means that 3 participants not infected with HIV reported positive result. This is called a false positive.
- In addition, 4.7% of participants (42/900) failed to obtain a results.

### REFERENCES

- WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices, Geneva: World Health Organization; 2019.
- Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
- WHO. TSS-1 Human Immunodeficiency Virus(HIV) rapid diagnostic tests for professional use and/or selftesting, Geneva: World Health Organization; 2016.

### SYMBOLS KEY

	In vitro diagnostic medical device		Consult instructions for use		Expiry date YYYY-MM-DD
	Contains sufficient for < n > tests		Date of manufacture YYYY-MM-DD		Keep dry
	Batch code		Temperature limit		Keep away from sunlight
	Manufacturer		Do not re-use		Product code/ Catalogue number
	Caution				

## Manufacturer information

**Guangzhou Wondfo Biotech Co., Ltd.**

**Add:** No.8 Lizhishan Road, Science City, Luogang District, 510663,Guangzhou, P.R. China  
**Tel:** +86-20-32296083 400-888-5268 (Toll Free)  
**Fax:** +86-20-32296063  
**E-mail:** global@wondfo.com.cn  
**Website:** www.wondfo.com.cn

Please contact the manufacturer or your local distributor if you have any questions related to the product.

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#### RÉFÉRENCES

1. WHO. TGS-5 Designing instruction for use for *in vitro* diagnostic medical devices. Geneva: World Health Organization ; 2019.
2. Medicines and Healthcare products Regulatory Agency. Guidance for notified bodies on the regulation of IVDs for self-testing. London, U.K.: MHRA, Competent Authority (UK); 2012.
3. WHO. TSS-1 Human Immunodeficiency Virus (HIV) rapid diagnostic tests for professional use and/or selftesting. Geneva: World Health Organization; 2016.

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