**THE UNITED REPUBLIC OF TANZANIA**

**MINISTRY OF HEALTH**

**TANZANIA MEDICINES AND MEDICAL DEVICES AUTHORITY**

**PUBLIC ASSESSMENT REPORT FOR CARESTART™ MALARIA PF/PAN(HRP2/PLDH) AG  
COMBO RDT**

**Version number 2.0, 09/09/2024**

**TMDA Headquarters, Plot No. 56/1, Block E, Kisasa B Centre, Swaswa Road, P. O. Box  
1253, Dodoma – Tanzania, Telephone: +255 (26) 2961989/2061990/+255 (22)  
2450512/2450751/2452108, Email: [info@tmda.go.tz](mailto:info@tmda.go.tz), Website: [www.tmda.go.tz](http://www.tmda.go.tz)**

**Toll free: 0800110084**

## 1. Introduction

*CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT* is a class C *in-vitro* diagnostic device belonging to the *Immunology* specialty category. It is approved in Tanzania as a *kit* for use in *adults, children, and elderly* by well-trained healthcare professionals.

### 1.1. Administrative Information

|                                                              |                                                                                                                                                                                                     |
|--------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Registration number                                          | TAN 24 MDR 0156                                                                                                                                                                                     |
| Brand Name (if relevant)                                     | CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT                                                                                                                                                  |
| Common name                                                  | Malaria Pf/Pan Ag                                                                                                                                                                                   |
| Class of the device and rule applied                         | Class C, Rule 3                                                                                                                                                                                     |
| GMDN code and term                                           | 52311: Multiple Plasmodium species antigen IVD, kit, rapid ICT, clinical                                                                                                                            |
| Name and complete address of the Market Authorization Holder | Access Bio, Inc. (Ethiopian Branch)<br>Addis Ababa, Yeka sub city, Woreda 9, Kebele 16,<br>House No-New, P.O. Box 22590, Ethiopia                                                                   |
| Name and address (es) of local responsible person (LRP).     | Reno Distributors Medical And Surgical Equipment,<br>P.O. Box 9202,<br>Swahili/Kipata Street Plot No 13, Block Number 46<br>House Number 37 Opposite To B&M Pharmacy<br>Dar es Salaam,<br>Tanzania. |

### 1.2. Assessment Procedure

The application for registration of *CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT* was submitted on 25/04/2024. The product underwent *abridged* assessment. Assessment was completed in two rounds of evaluation. *CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT* was registered on 29/08/2024.

## 2. Technical information

### 2.1. Intended use

The intended use of *CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT* as declared by the manufacturer and approved by TMDA is *for qualitative detection of malaria HRP2 (histidine-rich protein 2) of Plasmodium falciparum and pLDH*

(*plasmodium lactate dehydrogenase*) of *P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale* and *Plasmodium malariae* in human whole blood as an aid in the diagnosis of malaria infection. CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is approved for use in the laboratory or in point-of-care by well-trained healthcare professionals.

## 2.2. Device details and features

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is a single device with additional components or accessories. CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT has been registered as a kit which consists of cassettes (individually packaged in an aluminum pouch with a bag of desiccant), accessories including specimen transfer device, alcohol swab and lancet sufficient for testing, instructions for use (IFU) and buffer in a bottle.

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is a rapid test device. It is used for aid diagnosis of Malaria. CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT operates by detection of malaria HRP2 (histidine-rich protein 2) of *Plasmodium falciparum* and pLDH (*plasmodium lactate dehydrogenase*) of *P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale* and *Plasmodium malariae* in human whole blood. The RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN speci\_c to pLDH of the *Plasmodium* species (*P.falciparum*, *P.vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a Monoclonal antibody species to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles. The test out-put is qualitative.

The type of specimen used is whole blood and is collected by venous or capillary blood

CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT containing a test strip for lateral flow immunochromatographic assay encased in the cassette. The test strip is produced by cutting a component plate into strips. The plate is composed of a nitrocellulose membrane, a gold conjugate pad (pad dispensed with gold conjugate mixture - composed of monoclonal anti-pLDH antibody (6C9), monoclonal anti-HRP2 antibody (C1-13) and biotinylated monoclonal anti-pLDH antibody (19G7)), a filter pad and an absorbent pad attached on a plastic plate by adhesive material.

### Multi Kit

Multi kit is composed of CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT cassettes (individually packaged in an aluminum pouch with a bag of desiccant), accessories including specimen transfer device, alcohol swab and lancet sufficient for testing, instructions for use (IFU) and buffer in a bottle.

## Single Kit

Each single kit is consisted of a primary packaged CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT cassette and an accessory pack in a bag. 3.4.2 An accessory pack contains a specimen transfer device, alcohol swab, lancet, and buffer vial and instruction card packaged in a plastic bag. The content of accessory pack is subject to change upon the customer's request

### 2.3. Commercial presentation

There are 4 approved commercial presentations as follows:

| Consumables:<br>Configuration | Pack<br>size<br>(tests/kit) | Component                      |
|-------------------------------|-----------------------------|--------------------------------|
| Multi kit                     | 25                          | 1 × Assay buffer bottle (4 mL) |
|                               |                             | 25 × Specimen transfer device  |
|                               |                             | 25 × Lancet 2                  |
|                               |                             | 25 × Alcohol swab              |
|                               |                             | 1 × Instructions for Use       |
| Multi kit                     | 50                          | 1 × Assay buffer bottle (6mL)  |

|            |    |                                                                                                                                                                                          |
|------------|----|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|            |    | 50 × Specimen transfer device                                                                                                                                                            |
|            |    | 50 × Lancet 50 × Alcohol swab                                                                                                                                                            |
|            |    | 1 × Instructions for Use                                                                                                                                                                 |
| Single kit | 25 | Each single kit bag contains one of each of the following items: - Single use buffer vial - Specimen transfer device - Lancet - Alcohol swab - Instruction card 1 × Instructions for Use |
| Single kit | 40 | Each single kit bag contains one of each of the following items: - Single use buffer vial - Specimen transfer device - Lancet - Alcohol swab - Instruction card 1 × Instructions for use |

#### 2.4. Items required but not submitted

- a) Pair of disposable gloves
- b) Timer
- c) Sharps box
- d) Pencil or pen
- e) Sterile gauze or cotton

### 3. Storage instructions

#### 3.1.1. Shelf-life

The approved shelf-life is 30 months.

#### 3.1.2. Storage conditions

The recommended storage conditions is 1 to 40° C

#### 3.1.3. Shipping conditions

The recommended shipping conditions are *between 7.5°C and 35.5°C for 4 hours, and humidity between 45% and 66%.*

### 4. Manufacturing site audit

The manufacturer of the device is Access Bio, Inc. Ethiopian Branch 22590, [abe\\_bus@accessbio.net](mailto:abe_bus@accessbio.net) Addis Ababa Ethiopia.

The quality audit has not been conducted yet, but the application and payment have already been completed.

## 5. Performance Evaluation

### 5.1. Analytical Performance

The analytical performance characteristic of the device *was not submitted*.

### 5.2. Clinical Performance

Clinical performance was conducted at Access Bio, Inc. Ethiopian Branch 22590, [abe\\_bus@accessbio.net](mailto:abe_bus@accessbio.net), Addis Ababa Ethiopia. The sensitivity and specificity parameters were tested.

Based on results of the performance studies, it was concluded that the *test sensitivity are 98% and 96% for plasmodium falciparum and plasmodium vivax respectively. The specificity is 97.5% for both species.* The studies further concluded that CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is capable of consistently producing accurate and reliable test output.

## 6. Product label and instructions for use


The content of the primary and secondary labels is in line with TMDA labeling requirements in terms of content, layout and design. The label contains sufficient information for proper identification of the device and post marketing follow up of the product in the market.

The *instructions for use* include all the relevant information to ensure correct and safe use of the device by intended user.

### 6.1. Primary pack



## 6.2. Secondary pack

### 6.3 Instructions for use/Package insert

**CareStart™ Malaria Pf/PAN(HRP2/pLDH) Ag Combo RDT**

REF RMRU-02591 / RMRU-04091

#### TEST PROCEDURE

- Put on a new pair of gloves.
- Write the patient's name on the cassette.
- Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.
- 

Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.
- Wipe out the first drop of blood with sterile gauze or cotton.
- 

Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.

- 

Add 5µl of whole blood into "S" well.
- 

Add 3 drops (60µl) of buffer solution into "A" Well. Start a timer.
- Read result at 20 minutes.

#### INTERPRETATION OF THE TEST RESULT

**Negative**  
The presence of a line next to "C" indicates a negative result.

**Pf Positive**  
The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.

**PAN Positive**  
The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.

**Pf Positive or Mixed infection**  
The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale* or *P. malariae*).

**Invalid**  
The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.

For detailed information about the product and procedure, refer to the full version of instructions for use (Doc. No. IFU-RMRU71-EFSP) included in the RDT box.

IFU-RMRU91-E / Rev. D
May 21, 2018

## 7. Risk – Benefit Analysis

On basis of the data submitted, the current state of knowledge and compliance of the manufacturer to ISO 13485, the benefit of the product outweighs the risks associated with its use when used in accordance with the manufacturer instruction. *CareStart™ Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT* was recommended for registration.

## 8. Post-approval updates

### 8.1. Variation applications

| Reference number | Date submitted | Change requested | Recommendation | Granting date |
|------------------|----------------|------------------|----------------|---------------|
| NA               | NA             | NA               | NA             | NA            |

### 8.2. Feedback from pharmacovigilance, post marketing surveillance and enforcement activities

| Type of feedback          | Impact | Response |
|---------------------------|--------|----------|
| No recorded Adverse Event | NA     | NA       |

### 8.3. Re-registration applications

NA

## CHANGE HISTORY

| Version number | Date | Description of update | Section(s) Modified | Approval date |
|----------------|------|-----------------------|---------------------|---------------|
|                |      |                       |                     |               |

# CareStart™ Malaria Pf /PAN (HRP2/pLDH) Ag Combo RDT

## ENGLISH

**A RAPID DIAGNOSTIC TEST FOR THE DETECTION OF MALARIA HRP2 (Pf specific) AND pLDH (PAN specific) IN HUMAN BLOOD.**

## INTENDED USE:

For the rapid qualitative detection of malaria HRP2 (histidine-rich protein 2) of *Plasmodium falciparum* and pLDH (plasmodium lactate dehydrogenase) of *P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale* and *Plasmodium malariae* in human whole blood as an aid in the diagnosis of malaria infection.

## SUMMARY:

Malaria is a serious parasitic disease characterized by fever, chills and anemia. It is caused by a bite of infected Anopheles mosquitoes resulting in transmission of protozoan parasites to human. There are four major types of human malaria: *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae*. When infected, the parasites (called sporozoites) migrated to the liver where they mature and then are released to the blood stream of human infecting the red blood cells. Malaria infection occurs in more than 90 countries worldwide but is mostly prevalent in sub-Saharan Africa. It is estimated that there are over 300 million clinical cases and about 1 million malaria-caused deaths per year.

*CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the *Plasmodium* species (*P. falciparum*, *P.vivax*, *P. malariae* and/or *P. ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P. falciparum*. The conjugate pad is dispensed with antibodies adsorbed on gold particles.<sup>(1), 2), 4), 5)</sup>

*CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT is designed for the differentiated diagnosis of *P. falciparum*, *P.vivax*, *P. malariae* and *P. ovale* infection.

## WARNINGS AND PRECAUTIONS:

- Kits are for *in vitro* diagnostic and professional use only.
- Read the provided instructions for use before using the kit and follow the provided information when using the kit.
- Wash the hands thoroughly before and after using the kit.
- Wear protective gloves at all time while using the kits and dispose the gloves immediately after each testing.
- Change gloves and wash hands when contacted with potentially infectious materials.
- Do routine clean-up using an appropriate disinfectant.
- Do not eat or smoke while using the kit.
- Use the cassette and accessories (lancet and alcohol swab) immediately after opening the package.
- All provided materials are single-use. Do not re-use any of the materials.
- Do not use the materials from different lots.
- Do not use a cassette, if the packaging of cassette is compromised or the desiccant colour indicator is blue or black.
- Do not use the kit if the expiration date is past.
- The lancet is sterile. If a lancet cap is loose or damaged, do not use the lancet.
- Do not swallow the buffer solution.
- Do not use the alcohol swab if the package of the alcohol swab is damaged.
- Observe the storage condition (1~40°C) indicated on the packaging and RDT box.
- Do not freeze the kit but store the kit in a refrigerator. The refrigerated cassettes need to be warmed to the room temperature prior to use.
- Dispose wastes in accordance with the local regulations.
- Buffer bottle should be tightly closed after each use and stored in a cool area avoiding direct sun light. Discard the buffer solution and the bottle after the expiration date indicated on the bottle.

## PERFORMANCE CHARACTERISTICS:

The *CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT has been tested with positive and negative clinical samples confirmed by microscopic examination. The results are shown below:

| Specimen             | Positive | Negative | Sensitivity |
|----------------------|----------|----------|-------------|
| <i>P. falciparum</i> | 98       | 2        | 98%         |
| <i>P. vivax</i>      | 96       | 4        | 96%         |

| Specimen | Positive | Negative | Specificity |
|----------|----------|----------|-------------|
| Negative | 5        | 195      | 97.5%       |

95% CI, (97.05, 100) for *P. falciparum* positive result, (92.16, 99.84) for *P. vivax* positive result and (94.64, 99.36) for negative result (laboratory study)

## PRECISION:

Precision was evaluated by testing 10 replicates of *CareStart™* Malaria Pf/PAN (HRP2/pLDH) Ag Combo RDT with three specimens: a negative, a low positive and a strong positive specimen. The test results showed 100% confirmation of the expected results.

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Put on a new pair of gloves.
- 2

EN

Write the patient's name on the cassette.
- FR

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- 3

EN

Clean the area to be pierced using an alcohol swab. Let the alcohol dry completely before proceeding to the next step.
- 4

EN

Squeeze the end of a fingertip and pierce the cleaned area of the fingertip using a lancet provided. Discard the lancet in the sharps box.
- FR

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
- 6

EN

Collect the blood sample (5µl) using a provided specimen transfer device or a micropipette.
- 7

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- 8

EN

Add 2 drops (60µl) of buffer solution into "A" Well. Start a timer
- 9

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

Negative / Négatif / Negativo / Negativo

- EN

The presence of a line next to "C" indicates a negative result.
- FR

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Pf Positive / Pf Positif / Positivo para Pf / Pf Positivo

- EN

The presence of two lines (one line in the result window next to "C" and another in the result window next to "1") indicates a positive result for *P. falciparum*.
- FR

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- EN

The presence of two lines (one line in the result window next to "C" and another line in the result window next to "2") indicates a positive result for *P. falciparum*, *P. vivax*, *P. malariae* and/or *P. ovale*.
- FR

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- EN

The presence of three lines (3 lines in the result window next to "C", "2" and "1") indicates a positive result for *P. falciparum* or mixed infection of *P. falciparum* and other malaria (*P. vivax*, *P. ovale*, or *P. malariae*).
- FR

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- EN

The test is invalid when a line does not appear next to "C". If this occurs, the test should be repeated using a new cassette.
- FR

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LIMITATIONS AND INTERFERENCES:

- The following anticoagulants have been validated for use with this test: heparin, EDTA and citrate.
- This test is designed to detect HRP2 and pLDH antigens of Malaria *Plasmodium* species. Other clinically available tests are required if the obtained results are questionable. A definitive clinical diagnosis should not be made solely based on the result of this test, but should only be made by a qualified physician after all clinical and laboratory findings have been evaluated.
- A faint test line should be read as positive. A false-negative is possible due to a low parasite density.
- This test may still produce a positive result after successful anti-malarial treatment. Therefore, its use is not recommended for monitoring a response to anti-malarial treatment.<sup>3)</sup>
- If specimens cannot be tested immediately, they should be refrigerated at 2 - 8°C for up to 3 days.
- The test may produce a false positive result for a patient with acute schistosomiasis, a high level of rheumatoid factor or presence of human anti-mouse IgG antibody.<sup>6)</sup>
- No cross reactivity has been observed when the samples from a patient of chagas, dengue and leishmaniasis.
- Results are stable for 30 minutes after assay time of 20 minutes.
- The prozone effect may cause a false-negative result.

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REFERENCES / RÉFÉRENCES / REFERENCIAS / REFERÊNCIAS

- Valecha N., Eapen A., Usha Devi C., Ravindran J., Aggarwal A., and Subbarao SK. (2002). Field Evaluation of the ICT Malaria P.f./P.v. Immunochromatographic Test in India. Annals of Tropical Medicine & Parasitology. 96, 333-336.
- Iqbal J., Hira PR., Sher A., and AL-Enezi AA. (2001). Diagnosis of Imported Malaria by Plasmodium Lactate Dehydrogenase (pLDH) and Histidine-Rich Protein 2 (PfHRP-2)-based Immunocapture Assays. American Journal of Tropical Medicine and Hygiene. 64, 20-23.
- Tijtra E., Suprianto S., Dyer ME., Currie BJ. and Anstey NM. (2001). Detection of Histidine-Rich Protein 2 and Panmalarial ICT Malaria P.f./P.v. Test Antigens After Chloroquine Treatment of Uncomplicated falciparum Malaria Does Not Reliably Predict Outcome in Eastern Indonesia. American Journal of Tropical Medicine and Hygiene. 65, 593-598.
- Panton LJ., McPhie P., Malloy WL., Welles TE., Taylor DW. and Howard RJ. (1989). Purification and Partial Characterization of an Unusual Protein of Plasmodium falciparum: Histidine-Rich Protein II. Molecular and Biochemical Parasitology. 35,149-160.
- Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras. (1995). Plasmodium falciparum and Plasmodium vivax: Lactate Dehydrogenase Activity and its Application for in vitro Drug Susceptibility Assay. Experimental Parasitology. 80, 260-271.
- Malaria Rapid Diagnostic Test Performance. Results of WHO product testing of malaria RDTs: round 7 (2015-2016). Table A4.7: Phase 2 false-positive rate for *P. falciparum* in samples containing specific non-malarial infectious pathogens. 124.

**Access Bio Ethiopia**  
Yeka Sub-City, Woreda 9, Kebele 16, House No-New, Addis Ababa Ethiopia  
Tel.: +251 116675030/10  
Email: info@accessbio.net  
Web: www.accessbio.net